

RAPID-C+

Real-time Air Microbial Detector for Aseptic Filling Lines

Dr. Svetlana Kiseleva, CPO at Plair SA
safonina@plair.ch

Monday, April 7, 2025



ENVIRONMENTAL MONITORING IN STERILE MANUFACTURING

Environmental monitoring (EM) determines the quality of a controlled environment via microbial data collection. Data comes from samples of air, surfaces, and personnel in a clean space.

Production lines like RABS (left) or Isolators (right) are designed to keep the process into a grade A environment while the line is situated into a lower grade environment (grade B for RABS, grade C for Isolator).



Total Particle counting

Active air sampling

Settling plates



TRADITIONAL GROWTH BASED METHODS

Viable particle monitoring is traditionally performed with:

Active air samplers

Settle plates



Benefits


- Established method
- Compendial description
- Accepted by authorities
- Simple methods
- No validation needed anymore
- Current GMP / Pharmacopeial limits are in CFU

Challenges

- Labor intensive, human intervention creates risks of contamination
- Slow, requires growth of the microorganisms, identification may take days to weeks
- Only culturable microorganisms are detected
- Elaborate handling
- In certain cases, recovery rate is low

BFPC – OPPORTUNITY OF AN ALTERNATIVE METHOD

Viable particle monitoring in real-time using Laser-Induced Fluorescence (LIF)



Benefits

- In line monitoring
- Real-time detection
- Continuous monitoring

Initial challenges

- High rate of false positives
- Lack of equivalency with traditional methods
- Not possible to decontaminate instrument

*BFPC – bio-fluorescent particle counter

ENVIRONMENTAL MONITORING AS REFLECTED IN EU GMP ANNEX 1

**As part of the contamination control strategy (CCS),
environmental monitoring should now ensure:**

Continuous sampling (9.24)

Microbial identification (9.31)

Rapid and automated methods are encouraged to improve the established ones

“The adoption of suitable alternative monitoring systems such as rapid methods should be considered by manufacturers in order to expedite the detection of microbiological contamination issues and to reduce the risk to product. These rapid and automated microbial monitoring methods may be adopted after validation has demonstrated their equivalency or superiority to the established methods.” (9.28)

OUR FLAGSHIP PRODUCT RAPID-C+

A real-time air microbial detector for environmental monitoring in sterile production



Viable
Particle
Counts

Total
Particle
Counts


Continuous
active
sampling

All you need for grade A, outside of grade A

OPTICAL STAGE



Advanced
BFPC




Viable
Particle
Counts

(0.5-100 μm)

Total
Particle
Counts

In accordance
with ISO 21501-4

SAMPLING STAGE



No air concentrator, all particles monitored and impacted

C-CATCH ACTIVE SAMPLER

- C-Catch is a patented aseptic container to protect agar plate (AP)
- C-Catch is a certified active sampler (EN 17141:2020)
- Outside of grade A for quick and easy replacement of APs
- Direct compliance with CFUs
- Direct correlation with viable particles count possible



Sieve




AP ø90 mm



C-Catch

THE VERSATILITY OF RAPID-C+



TPC – Total particle count
Real-time and quantified


VPC – Viable particle count
Real-time and quantified

CFU – Colony Forming Unit
Continuous and quantified
(external incubation required)


All monitoring done on the same m³ of air

INTEGRATION INTO FILLING LINES

- Ready for integration into fill finish lines
- One exhaust for operation (internal pump)
- Second exhaust for decontamination (external pump)
- Active sampling qualified with complete tubing
- Easy access to replace AP



RAPID-C+ INTEGRATION INTO FILLING LINES



FEATURES

- Continuous viable and total particle monitoring
- Real-time and digitalized data

BENEFITS


- Faster batch release
- More sustainable production with less waste

COMMUNICATION INTERFACE


Rapid-C+ is a “probe” that can be interrogated via REST API (user, admin, service access levels)

Rapid-C+ REST API specification enables the coding of drivers in HMI / EM SW Platform


Plair Dashboard is a diagnostic tool with data/status-control/display



RAPID-C+ ANSWERS



THE BENEFITS OF USING RAPID-C+



Thank you !



Route de Saint-Julien 275,
1258 Perly-Certoux (Geneva), Switzerland
info@plair.ch / +41 (22) 552-3830

www.plair.ch



RAPID-C+
PLAIR

USER INTEREST GROUP ANNUAL MEETING

11 June 2025 | Geneva, Switzerland



Shaping the Future of BFPC and
Aseptic Filling: Join the Experts at
This Year's Must-Attend Event

